

Micropipette holder

Tool holder


Dimensions 18 x 6 x 7.5 mm
(Optimized for 55 mm micropipette)

Max. temperature 120° / ≈1h30 (autoclave)

Tubing

Dimensions and material Length: 2 m
Internal diameter: 1.3 mm
Material: Tygon®

External tube holder 80 x 31.5 x 1.5 mm



Included items:

- 1 x tube
- 1 x tube holder
- 2 x M6 screws
- 2 x rubber o-rings

Compatibility:

- Glass capillaries OD ≈ 1mm
- Third party injectors

